


## Using a microscope

Draw lines from the boxes to the diagram to show the parts of the microscope.



**Eyepiece Lens:** the lens at the top that you look through

**Stage:** the flat platform where you place your slides

### Method (What we did)

First, we \_\_\_\_\_ the slide onto the stage of the microscope.  
\_\_\_\_\_ we \_\_\_\_\_ the light to illuminate the onion cells. Next we \_\_\_\_\_ the knob so the second lens almost touched the slide. After this we \_\_\_\_\_ the eyepiece lens. \_\_\_\_\_ we \_\_\_\_\_ the knob \_\_\_\_\_ to bring the onion cells into focus.

*Use the words in the box below to **fill-in** the blanks in the text/paragraph/passage above.*

or

*Using the words in the box below **complete** the passage above. Use each word only once - **cross off** each word once it has been used.*

Finally	turned	upwards	looked down
put	switched on	Then	turned

*Write instructions on how to use a microscope in the correct sequence.*

Turn upwards	Put
Switch on	Look down
to illuminate	

Instruction ( What to do)	Correct order (Sequence)
A	
B	
C	
D	
E	